

My FibroScan
examination appointment:

Date: _____

Time: _____

Place: _____

Results obtained with
FibroScan will help
your physician manage
your liver health.

echosens

because liver health matters

Echosens North America
1 Research Drive, Suite 400B
Westborough, MA 01581
781-790-0845

FibroScan® Indications for Use: FibroScan is intended to provide 50Hz shear wave speed measurement for a wide range of tissue stiffness as well as 2.5 MHz transient elastography (TE) information (CAP: Controlled Attenuation Parameter) for intra-hepatic structure of the body. FibroScan is indicated for non-invasive assessment of liver fibrosis, liver stiffness, wave speed and estimates of stiffness, as well as determining a 3.5 MHz attenuation coefficient for the attenuation parameter (CAP: Controlled Attenuation Parameter). The shear wave speed of stiffness and CAP may be used as an initial diagnosis and monitoring of adult patients with liver disease, as part of the overall management of the liver. Shear wave speed of stiffness may be used as an initial or follow-up management of pediatric patients with liver disease.

© 2023 Echosens. All rights reserved. Echosens, FibroScan®, and all other marks are trademarks of Echosens in the United States and are registered in the U.S. and/or other countries.


What to expect from a FibroScan® examination

A rapid, painless, and
completely non-invasive
exam of your liver

FibroScan
by echosens

What is FibroScan®?

Examination with FibroScan®, also called transient elastography, is a non-invasive way to measure the stiffness of your liver. FibroScan® works by emitting a small pulse of energy, which may feel as a slight vibration on your skin. FibroScan® calculates the speed of this energy to give your healthcare provider an immediate measure of the stiffness of your liver. This stiffness measure can be an important part of understanding your overall liver health.



Is it painful?

No. The FibroScan® examination is painless, quick and easy. During measurement, you feel a slight vibration on the skin at the tip of the probe.



What will the FibroScan® examination be like?

You will lie on your back, with your right arm raised behind your head. Your healthcare provider will apply a water-based gel to your skin and then place the probe on the right side of your ribcage.


The healthcare provider will obtain a minimum of 10 quick measurements during your examination. These measurements will be used to establish an overall stiffness score.

FibroScan®
by echosens

The FibroScan® examination is covered
by most insurance plans.

Who are the best candidates for FibroScan®?

FibroScan® probes are adapted to work with most patients' body types and for all ages, from pediatric patients to adults.



How should I prepare for my FibroScan® exam?

It is important to have your FibroScan exam on an empty stomach. **Do not eat or drink anything for at least 4 hours prior** to your examination.

Wear comfortable clothes that will allow your healthcare provider to expose the right side of your rib cage.

To learn more, ask your physician or healthcare provider or visit:

echosens.com • myliverexam.com
patient.actonliverhealth.com

